


PMS-400A

800 MHz Gated Photon Counter / Multiscaler

- 2 Counter Channels
- 800 MHz Count Rate, 32 bit Resolution
- Direct Interfacing to most Detectors
- Multiscaler Mode: Up to 64k Time Channels, min. 250ns / Channel
- Gated Photon Counting: 1.5 ns min. Gate Pulse Width
- Event Recording Mode: Up to 32 k Events
- new** • 32 bit Accumulation Counter for ultra-fast Accumulation
- On-Board Discriminators, Timing and Control Logics
- new** • PCI Board with fast DMA (Bus Master),
- Software for Windows 98, NT, 2k and XP, Parallel Operation of Several Modules Supported



Becker & Hickl GmbH
Nahmitzer Damm 30
12277 Berlin
Tel. 030 / 787 56 32 <http://www.becker-hickl.com>
Fax. 030 / 787 57 34 email: info@becker-hickl.com

US Representative:
Boston Electronics Corp
91 Boylston Street, Brookline,
Massachusetts 02445 USA
Tel: (800) 347 5445 or (617) 566 3821
Fax: (617) 731 0935
www.boselec.com tcspc@boselec.com



Optical Transient Waveform Recording

The waveform of the light is measured with a resolution down to 250ns. Two signals can be recorded simultaneously. Applicable to luminescence decay of inorganic samples, phosphorescence, delayed fluorescence, chemoluminescence, LIDAR.

New: The PMS-400A provides a 32 bit accumulation counter which enables accumulation with virtually no dead time between sweeps (< 100 ns).

Recording of Luminescence Spectra

The luminescence and the excitation light are recorded simultaneously. Corrected excitation spectra are obtained by calculating B/A.

Single Molecule Detection

Recording of photon bursts. If the count rate inside a programmed time interval exceeds a programmed value, the number of photons and the time of the event is stored.

Gated Detection

The gate is opened during the laser pulse only. Events outside the laser pulses are suppressed. Exceptionally low background count rate.

Gating off Scattering Pulses

The gate is closed during the laser pulses. Scattered photons during the laser pulses are suppressed, the luminescence photons outside the laser pulses are recorded.

Specification (Typical Values)

Counter Channels	2
Count Rate (Input Amplitude 50mV, peak-peak)	800 MHz
min. Count Pulse Width	800 ps
min. Gate Width (Input Amplitude 200mV, peak-peak)	1 ns
min Trigger Pulse Width	1 ns
Discriminator Threshold (Count Inputs)	-1 V to +1 V in steps of 4 mV
Discriminator Threshold (Gate Inputs)	-2 V to +2 V in steps of 16 mV
Discriminator Threshold (Trigger Input)	-2 V to +2 V in steps of 16 mV
Input Connectors	MCX, 50 Ω
Counter Width	32 bit
Accumulation Counter	32 bit
Dead time between sweeps	< 100ns
No. of Time Bins	64 k for each counter channel
Time / Bin	250 ns to 100 000 s
Hardware Environment	Pentium PC
Software Environment	Windows 95, 98, 2000 or NT
Dimensions	180 mm x 108 mm x 15 mm

